


# EGIS Colorectal Stent

EGIS Colorectal Stent fits naturally into the lower digestive tract and allows patients to defecate normally. The proximal, middle and distal GOLD MARKERS increase visibility in fluoroscopy.



Double Covered Type

EGIS Colorectal Stent is used to open strictures in the large intestine. The Stent is designed with outstanding conformability to fit naturally into the large intestine.



## Specification

Delivery System 10Fr, 80/150/220cm | 12Fr, 80/150cm  
 Ordering No. First two\*\* : Diameter(mm) / Last two\*\* : Length(cm)

TYPE	Ordering No.		Size	
	OTW	TTS	Diameter	Length
Single Bare	CSB****08	CSB****22	18, 20, 22, 24, 26, 28, 30	6, 8, 10, 12
	CSB****15			
Double Bare	CDB****08	CDB****22		
	CDB****15			
Single Covered	CSC****08	CSC****22		
	CSC****15			
Double Covered	CDC****08	CDC****22		
	CDC****15			