


EGIS Pyloric Stent

EGIS Pyloric Stent is designed to fit naturally into the gastric outlet and duodenum. The proximal, middle and distal GOLD MARKERS increase visibility in fluoroscopy.



Pyloric Double Bare Type

EGIS Pyloric Stent is used for patients with a stricture or fistula on the gastric outlet and duodenum. The optimal radial force of the EGIS Pyloric stent will allow patients to achieve normal digestion. The stent is designed with excellent conformability to fit naturally into the gastric outlet and duodenum



Specification

Delivery System 10Fr, 120/180cm

Ordering No. First two** : Diameter(mm) / Last two** : Length(cm)

TYPE	Ordering No.		Size	
	OTW	TTS	Diameter	Length
Single Bare	PSB****12	PSB****18	18, 20, 22, 24	6, 8, 10, 12
Double Bare	PDB****12	PDB****18		
Single Covered	PSC****12	PSC****18		
Double Covered	PDC****12	PDC****18		