



# Reboa Ballon Kit™

[www.reboamedical.com](http://www.reboamedical.com)

World's first complete REBOA kit  
In/Out of hospital  
Trauma  
PPH

Kit contains; Reboa Balloon, drape, introducer, guidewire, and more. For more information, contact your local distributor or Reboa Medical directly.


Product code	Balloon (Ø/L)	Fill volume	Introducer Size (included in kit)
RBK15305006	15/30 (mm)	8,0 mL	6F
RBK20305007	20/30 (mm)	15,0 mL	7F

\*both kits are delivered with a 23 cm long introducer

The REBOA balloon is inserted using standard Seldinger technique.



Approved by CE, FDA, and Health Canada.



**“PROBLEMS CAN BE COMPLICATED.  
SOLUTIONS CAN NOT”**