

SB Hood


Hard type series (For Upper GI)



Soft type series (For Upper and Lower GI)




Short type




Bottom projection type



Straight type




Undercut type



Top projection type




Features



SB Hood (Hard) short type

Short type provides wider endoscopic view and secure procedures.



SB Hood (Hard) Bottom projection type

By pressing Bottom projection type towards the mucosal layer, the dissection area can be seen in the center of endoscopic view.




SB Hood (Hard) Top projection type

Top projection type slides into the submucosal layer and produces counter traction for secure dissection.




Features

*SB Hood (Soft)
Straight type*




*SB Hood (Soft)
Undercut type*



3.5mm

Undercut type

(more capable to slide under tissue)



- Very soft and can be used for ESD, especially colon.
- Clear material to obtain good endoscopic view.
- Uniquely designed side hole for ventilation.
- Inclined tiptop to approach under tissue easily.
- 6mm length & oval shaped front hole.
(makes SB knife Jr & Short type more effective)

Model number	Type	Length of Protrusion from Endoscope	Maximum O.D.	Suitable O.D. of Endoscope	Colour of Joint	Packaging
MD-47910	Hard Top Projection(TP)	12mm	13.3mm	9.9mm	Clear	12pieces / case
MD-47920	Hard Bottom Projection(BP)	12mm	13.3mm	9.9mm	Blue	12pieces / case
MD-47930	Hard Short	4mm	13.3mm	9.9mm	Clear	12pieces / case
MD-47940	Hard Set of 3 hard types	-	13.3mm	9.9mm	-	12pieces / case
MD-47950	Soft Straight	6mm	13.0mm	9.9mm - 10.8mm	-	5 pieces / case
MD-47951	Soft Under cut	Top 6mm Bottom 3.5mm	13.0mm	9.9mm - 10.8mm	-	5 pieces / case
MD-47960	Soft Straight	6mm	14.8mm	11.7mm -12.3mm	-	5 pieces / case
MD-47961	Soft Under cut	Top 6mm Bottom 3.5mm	14.8mm	11.7mm -12.3mm	-	5 pieces / case
MD-47970	Soft Straight	6mm	15.9mm	12.8mm -13.2mm	-	5 pieces / case
MD-47971	Soft Under cut	Top 6mm Bottom 3.5mm	15.9mm	12.8mm -13.2mm	-	5 pieces / case



SUMITOMO BAKELITE CO.,LTD.

<http://www.sumibe.co.jp>

Global Sales Dept.

Medical Products Business Div.

2-5-8, Higashi-Shinagawa, Shinagawa-ku, Tokyo, 140-0002, Japan

Tel: 81-3-5462-4824

Fax: 81-3-5462-4894