



VenusP-Valve™ System – Model Selection Guide

Dimension of VenusP-Valve™ when fully expanded



Model	Specification	Diameter [mm]			D [mm]
		A	B	C	
L28P	P28-25	38.0	28.0	38.0	25.0
	P28-30				30.0
L30P	P30-25	40.0	30.0	40.0	25.0
	P30-30				30.0
L32P	P32-25	42.0	32.0	42.0	25.0
	P32-30				30.0
L34P	P34-25	44.0	34.0	44.0	25.0
	P34-30				30.0
L36P	P36-25	46.0	36.0	46.0	25.0
	P36-30				30.0



Dimension of Delivery Catheter System

Model	Specification	Delivery Catheter System [mm]	
		Capsule OD	Loading Length
22Fr	DS22-P77	7.33	77
	DS22-P70		70
24Fr	DS24-P87	8.00	87
	DS24-P81		81
	DS24-P75		75

Device components matching table. VenusP-Valve™ with matching Delivery Catheter System

Percutaneous Pulmonary Valve		Delivery Catheter System	
Model	Specification	Model	Specification
L28P	P28-25	22FR	DS22-P70
	P28-30		DS22-P77
L30P	P30-25		DS22-P70
	P30-30		DS22-P77
L32P	P32-25	24FR	DS24-P75
	P32-30		DS24-P81
L34P	P34-25		DS24-P81
	P34-30		DS24-P87
L36P	P36-25		DS24-P81
	P36-30		DS24-P87

Product not intended to enter the market in Germany, Ireland or France.
 VENUSMEDTECH, VenusP-Valve, the stylized QI Logo are all trademarks of Venus Medtech Inc.
 © 2022 Venus Medtech Inc. All rights reserved.