


For CRYO-S ELECTRIC II

CRYOPROBES CATALOGUE


▪ ophthalmology



how to identify the probe?

order/ref. no.

0 2 3004 0



In care for highest quality of offered products and services, Metrum Cryoflex reserves the right to modify products and introduce changes in descriptions.

Ophthalmology

Cryosurgery in ophthalmology has a rich history and continues to be an important supplement in the treatment of ophthalmic pathology. The use of cryotherapy in ophthalmology has helped advance maturing fields such as cataract extraction, while in other instances revolutionized patient care including, trichiasis, retinopathy of prematurity, ingrown eyelash and ocular surface malignancies.

Standard accessories	Qty
2-button footswitch with a cable	1
Cylinder-unit pneumatic hose	2
Flat wrench (N ₂ O and CO ₂ connector)	1
Operating manual	1
Electric power cord	1
N ₂ O and CO ₂ cylinder connectors	2
Spare seals and seal cones	set
Spare fuse	2
Unit packaging	1
Socket caps (gas supply, cryoprobe)	2
Screw-in cryoprobe pantograph	1
Socket caps (gas supply, cryoprobe)	2
Screw-in cryoprobe pantograph	1
Moisture filter	set
Hose for detachable cryoprobe	1











CE 2274




Cryo-S Electric II

order/ref. no. **CSEII**

Accessories and spare parts	Ref. / Order No.
Trolley	TRCSE
Reusable cryoprobes	see catalog
Hose for detachable cryoprobe	HDPCSE
Sterilization cryoprobe plug (black)	CPCSE
Sterilization hose plugs set (black)	HPCSE
Moisture filter insert	FI
Spare seals set	SSCSE
Gas – CO ₂ – connector	COC
Gas – N ₂ O – connector	NOC 9
Cylinder – unit pneumatic hose (50cm)	HC5
Cylinder – unit pneumatic hose (100cm)	HC10



cylindrical cryoprobe $\varnothing 4/4$ **type: OP-H015S/RF****reusable, detachable****order/ref. no. 0 2 3004 1****reusable, fixed****order/ref. no. 0 2 3004 0**cylindrical cryoprobe $\varnothing 4/6$ **type: OP-H015L/RF****reusable, detachable****order/ref. no. 0 2 3006 1****reusable, fixed****order/ref. no. 0 2 3006 0**micro cryoprobe $\varnothing 1.0\text{mm}$ **type: OP-L008/RF****reusable, detachable****order/ref. no. 0 2 4010 1****reusable, fixed****order/ref. no. 0 2 4010 0**micro cryoneedle $\varnothing 0.7$ (for ingrown eyelashes)**type: OP-N008/RF****reusable, detachable****order/ref. no. 0 2 5008 1****reusable, fixed****order/ref. no. 0 2 5008 0**



micro cryoprobe for retina $\varnothing 1.6\text{mm}$



type: OP-R009/RF

reusable, detachable

order/ref. no. **0 2 6016 1**

reusable, fixed

order/ref. no. **0 2 6016 0**



micro cryoprobe for retina $\varnothing 3.0\text{mm}$


type: OP-R010S/RF

reusable, detachable

order/ref. no. **0 2 6130 1**

reusable, fixed

order/ref. no. **0 2 6130 0**



micro cryoprobe for retina $\varnothing 3.0\text{mm}$


type: OP-R010L/RF

reusable, detachable

order/ref. no. **0 2 6131 1**

reusable, fixed

order/ref. no. **0 2 6131 0**



short and angular cryoprobe for retina $\varnothing 3.0\text{mm}$



type: OP-R011S/RF

reusable, detachable

order/ref. no. **0 2 6230 1**

reusable, fixed

order/ref. no. **0 2 6230 0**



long and angular cryoprobe for retina $\varnothing 3.0\text{mm}$


type: OP-R011L/RF

reusable, detachable

order/ref. no. **0 2 6231 1**

reusable, fixed

order/ref. no. **0 2 6231 0**



cryoprobe for cataract $\varnothing 1.5\text{mm}$


type: OP-C013/RF

reusable, detachable

order/ref. no. **0 2 1013 1**

reusable, fixed

order/ref. no. **0 2 1013 0**



cryoprobe for glaucoma $\varnothing 3.0\text{mm}$




type: OP-G014/RF

reusable, detachable

order/ref. no. **0 2 2014 1**

reusable, fixed

order/ref. no. **0 2 2014 0**



flat cryoprobe for retina $2 \times 4\text{mm}$


type: OP-R102/RF

reusable, detachable

order/ref. no. **0 2 6440 1**

reusable, fixed

order/ref. no. **0 2 6440 0**



cryoprobe for the vitreous body $\varnothing 1.0\text{mm}$

type: OP-V008/RF

reusable, detachable

order/ref. no. **0 2 7010 1**

reusable, fixed

order/ref. no. **0 2 7010 0**



METRUM CRYOFLEX

Manufacture

Zielna 29 Street
05-082 Blizne Łaszczyńskiego
POLAND, EU

Headoffice

Kolejowa 16A Street
05-092 Łomianki
POLAND, EU

Tel: +48 22 33 13 750
+48 22 33 13830
Fax: +48 22 33 13 766

export@metrum.com.pl
www.metrum.com.pl