

Galt|STICK®



Embedded radiopaque
marker band on all sheaths



Echogenic,
marked Chiba needle



Echogenic,
marked Trocar needle

The GaltStick® Introducer System is a 6 French Coaxial sheath assembly with stiffening stylet. The GaltStick® allows for the introduction of a .018" guidewire and .038" working guidewire in percutaneous procedures.

- Embedded radiopaque marker band on all sheaths for enhanced visibility
- Smooth dilator/sheath transition
- Variety of kit offerings
- Hydrophilic or silicone coated sheath configurations
- 6Fr outer sheath, 4Fr dilator, 20cm usable length
- Echogenic, marked Trocar needle or Echogenic, marked Chiba needle



Silicone Coated GaltStick® Introducer Set with Radiopaque Band

Model Number	Hydrophilic	.018" x 60cm SS/ Gold guidewire	.018" x 60cm NT/ Gold guidewire	.038" x 145 PTFE 3mm "J" Straight guidewire	Echo Needle	Quantity
KIT-030-04		X		X	Trocar	5
KIT-030-05			X	X	Trocar	5
KIT-030-06		X			Trocar	10
KIT-030-03			X		Trocar	10

Hydrophilic Coated GaltStick® Introducer Set with Radiopaque Band

Model Number	Hydrophilic	.018" x 60cm SS/ Gold guidewire	.018" x 60cm NT/ Gold guidewire	.038" x 145 PTFE 3mm "J" Straight guidewire	Echo Needle	Quantity
KIT-036-04	X	X		X	Trocar	5
KIT-036-05	X		X	X	Trocar	5
KIT-036-06	X	X			Trocar	10
KIT-036-03	X		X		Trocar	10
KIT-037-03	X		X		Chiba	10
KIT-037-06	X	X			Chiba	10

Accessories

Model Number	Description	Quantity
NDL-101-00	21ga x 15cm Trocar, Echogenic	5
NDL-109-00	22ga x 15cm Chiba, Echogenic	5
NDL-109-01	22ga x 20cm Chiba, Echogenic	5

