

Tokai Microcatheter

Carnelian NT

NON-TAPER

Tri-axial Microcatheter System
with CARNELIAN 2.7

Tri-axial Microcatheter System

CE
1639

Carnelian NT

Carnelian 2.7

- Boost Support
- Distal Clarity
- Ease of Access
- Prevent Glue Migration
- Avoid Catheter Re-selection

Catalogue No.	Proximal (F)	Distal (F)	Length (cm)	Inner Diameter Proximal/Distal (inch)	Tip	Hydrophilic Coating (cm)	Guidewire (inch)	Guide Catheter (inch)
---------------	--------------	------------	-------------	---------------------------------------	-----	--------------------------	------------------	-----------------------

CARNELIAN NT - NON TAPER

MXN135X-E	1.9	1.9	135	0.017/0.017	Straight	80	0.014	0.027
MXN155X-E	1.9	1.9	155	0.017/0.017	Straight	110	0.014	0.027
MXNS135X-E	1.8	1.6	135	0.0155/0.0155	Straight	80	0.014	0.027
MXNS155X-E	1.8	1.6	155	0.0155/0.0155	Straight	110	0.014	0.027

CARNELIAN 2.7 - HIGH FLOW

MCH105-E	2.9	2.7	105	0.027/0.027	Straight	60	0.025	0.041
MCH125-E	2.9	2.7	125	0.027/0.027	Straight	60	0.025	0.041