



# EGIS Flower Stent




01



## EGIS Flower Stent

- Effective drainage of bile juice by side grooves
- Prevent cholecystitis and pancreatitis
- e-PTFE fully covered and retrieval stent


## Upgraded Special Features



- Innovative design with 5~6 grooves, facilitating bile drainage



- Fully covered stent with e-PTFE, designed for use in malignant & benign strictures



- Sufficient radial force to push out the stricture



- Low axial force, ensuring excellent flexibility & conformability



- For easy removal, a triangle-shaped lasso is attached



- To obtain anti-migration function, flared ends and grooves are made in the entire body

Delivery System 8Fr, 50/180cm

Ordering NO. First two\*\* : Diameter(mm) / Last two\*\* : Length(cm)

| Ordering No.  |              |              | Size     |                             |
|---------------|--------------|--------------|----------|-----------------------------|
| Type          | OTW          | TTS          | Diameter | Length                      |
| Bare          | BFBB****0805 | BFBB****0818 | 8, 10    | 4, 5, 6, 7,<br>8, 9, 10, 12 |
| Partial       | BFPB****0805 | BFPB****0818 |          |                             |
| Fully Covered | BFCB****0805 | BFCB****0818 |          |                             |



## EGIS Flower Stent

| Stent Type             | OTW(Over The Wire) |              |            |                 |            | TTS(Through The Scope) |              |            |                 |            |
|------------------------|--------------------|--------------|------------|-----------------|------------|------------------------|--------------|------------|-----------------|------------|
|                        | Ordering NO.       | Stent        |            | Delivery System |            | Ordering No.           | Stent        |            | Delivery System |            |
|                        |                    | Diameter(mm) | Length(cm) | Profile(fr.)    | Length(mm) |                        | Diameter(mm) | Length(cm) | Profile(fr.)    | Length(mm) |
| Flower Fully Covered   | BFCB08040805       | 8            | 4          | 8               | 500        | BFCB08040818           | 8            | 4          | 8               | 1800       |
|                        | BFCB08050805       |              | 5          |                 |            | BFCB08050818           |              | 5          |                 |            |
|                        | BFCB08060805       |              | 6          |                 |            | BFCB08060818           |              | 6          |                 |            |
|                        | BFCB08070805       |              | 7          |                 |            | BFCB08070818           |              | 7          |                 |            |
|                        | BFCB08080805       |              | 8          |                 |            | BFCB08080818           |              | 8          |                 |            |
|                        | BFCB08090805       |              | 9          |                 |            | BFCB08090818           |              | 9          |                 |            |
|                        | BFCB08100805       |              | 10         |                 |            | BFCB08100818           |              | 10         |                 |            |
|                        | BFCB08120805       |              | 12         |                 |            | BFCB08120818           |              | 12         |                 |            |
|                        | BFCB10040805       | 10           | 4          | 8               | 500        | BFCB10040818           | 10           | 4          | 8               | 1800       |
|                        | BFCB10050805       |              | 5          |                 |            | BFCB10050818           |              | 5          |                 |            |
|                        | BFCB10060805       |              | 6          |                 |            | BFCB10060818           |              | 6          |                 |            |
|                        | BFCB10070805       |              | 7          |                 |            | BFCB10070818           |              | 7          |                 |            |
|                        | BFCB10080805       |              | 8          |                 |            | BFCB10080818           |              | 8          |                 |            |
|                        | BFCB10090805       |              | 9          |                 |            | BFCB10090818           |              | 9          |                 |            |
|                        | BFCB10100805       |              | 10         |                 |            | BFCB10100818           |              | 10         |                 |            |
|                        | BFCB10120805       |              | 12         |                 |            | BFCB10120818           |              | 12         |                 |            |
| Flower Partial Covered | BFPB08040805       | 8            | 4          | 8               | 500        | BFPB08040818           | 8            | 4          | 8               | 1800       |
|                        | BFPB08050805       |              | 5          |                 |            | BFPB08050818           |              | 5          |                 |            |
|                        | BFPB08060805       |              | 6          |                 |            | BFPB08060818           |              | 6          |                 |            |
|                        | BFPB08070805       |              | 7          |                 |            | BFPB08070818           |              | 7          |                 |            |
|                        | BFPB08080805       |              | 8          |                 |            | BFPB08080818           |              | 8          |                 |            |
|                        | BFPB08090805       |              | 9          |                 |            | BFPB08090818           |              | 9          |                 |            |
|                        | BFPB08100805       |              | 10         |                 |            | BFPB08100818           |              | 10         |                 |            |
|                        | BFPB08120805       |              | 12         |                 |            | BFPB08120818           |              | 12         |                 |            |
|                        | BFPB10040805       | 10           | 4          | 8               | 500        | BFPB10040818           | 10           | 4          | 8               | 1800       |
|                        | BFPB10050805       |              | 5          |                 |            | BFPB10050818           |              | 5          |                 |            |
|                        | BFPB10060805       |              | 6          |                 |            | BFPB10060818           |              | 6          |                 |            |
|                        | BFPB10070805       |              | 7          |                 |            | BFPB10070818           |              | 7          |                 |            |
|                        | BFPB10080805       |              | 8          |                 |            | BFPB10080818           |              | 8          |                 |            |
|                        | BFPB10090805       |              | 9          |                 |            | BFPB10090818           |              | 9          |                 |            |
|                        | BFPB10100805       |              | 10         |                 |            | BFPB10100818           |              | 10         |                 |            |
|                        | BFPB10120805       |              | 12         |                 |            | BFPB10120818           |              | 12         |                 |            |
| Flower Bare            | BFBB08040805       | 8            | 4          | 8               | 500        | BFBB08040818           | 8            | 4          | 8               | 1800       |
|                        | BFBB08050805       |              | 5          |                 |            | BFBB08050818           |              | 5          |                 |            |
|                        | BFBB08060805       |              | 6          |                 |            | BFBB08060818           |              | 6          |                 |            |
|                        | BFBB08070805       |              | 7          |                 |            | BFBB08070818           |              | 7          |                 |            |
|                        | BFBB08080805       |              | 8          |                 |            | BFBB08080818           |              | 8          |                 |            |
|                        | BFBB08090805       |              | 9          |                 |            | BFBB08090818           |              | 9          |                 |            |
|                        | BFBB08100805       |              | 10         |                 |            | BFBB08100818           |              | 10         |                 |            |
|                        | BFBB08120805       |              | 12         |                 |            | BFBB08120818           |              | 12         |                 |            |
|                        | BFBB10040805       | 10           | 4          | 8               | 500        | BFBB10040818           | 10           | 4          | 8               | 1800       |
|                        | BFBB10050805       |              | 5          |                 |            | BFBB10050818           |              | 5          |                 |            |
|                        | BFBB10060805       |              | 6          |                 |            | BFBB10060818           |              | 6          |                 |            |
|                        | BFBB10070805       |              | 7          |                 |            | BFBB10070818           |              | 7          |                 |            |
|                        | BFBB10080805       |              | 8          |                 |            | BFBB10080818           |              | 8          |                 |            |
|                        | BFBB10090805       |              | 9          |                 |            | BFBB10090818           |              | 9          |                 |            |
|                        | BFBB10100805       |              | 10         |                 |            | BFBB10100818           |              | 10         |                 |            |
|                        | BFBB10120805       |              | 12         |                 |            | BFBB10120818           |              | 12         |                 |            |