

## 01 EGIS Biliary Fully Covered Stent



- ▣ e-PTFE fully covered stent with a lasso for both benign and malignant diseases
- ▣ Excellent flexibility and conformability with low axial force
- ▣ Sufficient radial force to push the strictures out



## 02 Specification

Delivery System : 8Fr, 50/180cm

Ordering NO. First two\*\* : Diameter(mm) / Last two\*\* : Length(cm)

| Ordering No.                |              |              | Size     |                             |
|-----------------------------|--------------|--------------|----------|-----------------------------|
| Type                        | OTW          | TTS          | Diameter | Length                      |
| Biliary Fully Covered Stent | BEPB****0805 | BEPB****0818 | 8, 10    | 4, 5, 6, 7,<br>8, 9, 10, 12 |

| Stent Type                     | OTW(Over The Wire) |              |            |                 |            |
|--------------------------------|--------------------|--------------|------------|-----------------|------------|
|                                | Ordering NO.       | Stent        |            | Delivery System |            |
|                                |                    | Diameter(mm) | Length(cm) | Profile(fr.)    | Length(mm) |
| Biliary Fully Covered (e-PTFE) | BEPB08040805       | 8            | 4          | 8               | 500        |
|                                | BEPB08050805       |              | 5          |                 |            |
|                                | BEPB08060805       |              | 6          |                 |            |
|                                | BEPB08070805       |              | 7          |                 |            |
|                                | BEPB08080805       |              | 8          |                 |            |
|                                | BEPB08090805       |              | 9          |                 |            |
|                                | BEPB08100805       |              | 10         |                 |            |
|                                | BEPB08120805       |              | 12         |                 |            |
|                                | BEPB10040805       | 10           | 4          | 8               | 500        |
|                                | BEPB10050805       |              | 5          |                 |            |
|                                | BEPB10060805       |              | 6          |                 |            |
|                                | BEPB10070805       |              | 7          |                 |            |
|                                | BEPB10080805       |              | 8          |                 |            |
|                                | BEPB10090805       |              | 9          |                 |            |
|                                | BEPB10100805       |              | 10         |                 |            |
|                                | BEPB10120805       |              | 12         |                 |            |

| Stent Type                     | TTS(Through The Scope) |              |            |                 |            |
|--------------------------------|------------------------|--------------|------------|-----------------|------------|
|                                | Ordering No.           | Stent        |            | Delivery System |            |
|                                |                        | Diameter(mm) | Length(cm) | Profile(fr.)    | Length(mm) |
| Biliary Fully Covered (e-PTFE) | BEPB08040818           | 8            | 4          | 8               | 1800       |
|                                | BEPB08050818           |              | 5          |                 |            |
|                                | BEPB08060818           |              | 6          |                 |            |
|                                | BEPB08070818           |              | 7          |                 |            |
|                                | BEPB08080818           |              | 8          |                 |            |
|                                | BEPB08090818           |              | 9          |                 |            |
|                                | BEPB08100818           |              | 10         |                 |            |
|                                | BEPB08120818           |              | 12         |                 |            |
|                                | BEPB10040818           | 10           | 4          | 8               | 1800       |
|                                | BEPB10050818           |              | 5          |                 |            |
|                                | BEPB10060818           |              | 6          |                 |            |
|                                | BEPB10070818           |              | 7          |                 |            |
|                                | BEPB10080818           |              | 8          |                 |            |
|                                | BEPB10090818           |              | 9          |                 |            |
|                                | BEPB10100818           |              | 10         |                 |            |
|                                | BEPB10120818           |              | 12         |                 |            |