




EGIS Biliary Flared Fully Covered (Valve) Stent



01 EGIS Flared fully covered stent




- Large Flare with 6mm difference from main body as an Anti-migration system (8mm stent with 14mm flare; 10mm stent with 16mm flare)
- Fully Covered type with Removable Lasso
- e-PTFE Cover providing outstanding biocompatibility
- Excellent flexibility and conformability due to the lowest possible axial force

02 EGIS Flared fully covered valve stent



- Large Flare with 6mm difference from main body as an Anti-migration system (8mm stent with 14mm flare; 10mm stent with 16mm flare)
- Fully Covered type
- e-PTFE Cover providing outstanding biocompatibility
- Excellent flexibility and conformability due to the lowest possible axial force
- Anti-Reflux Valve reduces the risk of post-ERCP cholangitis
- Upgraded Version from EGIS M-Valve stent; Fully Covered, Larger Flare, Stronger Radial Force, Lower Axial Force





EGIS Biliary Flared Fully Covered (Valve) Stent

Stent Type	Ordering NO.	Stent Type	Ordering NO.	Stent		Delivery System		
				Diameter(mm)	Length(cm)	Diameter(mm)	Usable Length(mm)	
flared fully covered stent	FFPC08040	flared fully covered valve stent	FVPC08040	14/8	4	2.7mm (8Fr)	500	
	FFPC08050		FVPC08050		5			
	FFPC08060		FVPC08060		6			
	FFPC08070		FVPC08070		7			
	FFPC08080		FVPC08080		8			
	FFPC08090		FVPC08090		9			
	FFPC08100		FVPC08100		10			
	FFPC08120		FVPC08120		12			
	FFPC10040		FVPC10040		16/10			4
	FFPC10050		FVPC10050					5
	FFPC10060		FVPC10060					6
	FFPC10070		FVPC10070					7
	FFPC10080		FVPC10080	8				
	FFPC10090		FVPC10090	9				
	FFPC10100		FVPC10100	10				
	FFPC10120		FVPC10120	12				
	FFEC08040		FVEC08040	14/8			4	
	FFEC08050		FVEC08050				5	
	FFEC08060		FVEC08060				6	
	FFEC08070		FVEC08070				7	
	FFEC08080		FVEC08080		8			
	FFEC08090		FVEC08090		9			
	FFEC08100		FVEC08100		10			
	FFEC08120		FVEC08120		12			
	FFEC10040		FVEC10040		16/10		4	
	FFEC10050		FVEC10050				5	
	FFEC10060		FVEC10060				6	
	FFEC10070		FVEC10070				7	
	FFEC10080		FVEC10080	8				
	FFEC10090		FVEC10090	9				
	FFEC10100		FVEC10100	10				
	FFEC10120		FVEC10120	12				

* Flared fully covered stent **Lasso : 20mm**