

TECHNICAL DETAILS

INTRODUCER SHEATH	
Material – outer	Stainless steel reinforced polymer
Material – inner	PTFE liner
Sheath coating	Hydrophobic, distal 15 cm
Sheath tip	Smooth taper with radiopaque platinum iridium marker band
French size	4, 5, 6 F
Length	45, 60, 90 cm
Sidearm/hemostatic valve	Removable, with color coded 3-way-stopcock


DILATOR	
Recommended guidewire	0.035"
Tip	Smooth taper, radiopaque
Shape	Straight

ORDERING INFORMATION

ARTICLE NR	DIAMETER (FRENCH)	LENGTH (CM)	SHAPE	UNITS/ CASE
FG-04045	4 F	45 cm	Straight	5
FG-04060	4 F	60 cm	Straight	5
FG-04090	4 F	90 cm	Straight	5
FG-05045	5 F	45 cm	Straight	5
FG-05060	5 F	60 cm	Straight	5
FG-05090	5 F	90 cm	Straight	5
FG-06045	6 F	45 cm	Straight	5
FG-06060	6 F	60 cm	Straight	5
FG-06090	6 F	90 cm	Straight	5

All parts suitable for guide wires up to 0.035 inch. Additional sizes to become available.

Manufactured by:
Contract Medical International GmbH
Lauensteiner Str. 37
01277 Dresden
Germany



KINK-RESISTANT
INTRODUCER SHEATH
Supporting complex access

SUPPORTING COMPLEX ACCESS

From easy insertion to smooth navigation to superb visibility:
DuraSheath is designed to combine optimal performance with excellent pushability and trackability. While maintaining the inner diameter, DuraSheath allows for smooth device passage and has an exceptional kink-radius and non-ovaling effect. The detachable hemostatic valve provides flexibility and versatility, giving you the ability to focus on the intervention.

