


# D'VILL™

Long-length, kink-resistant hemostasis  
introducer for introduction of large  
diameter interventional devices

Hemostasis valve  
designed to minimize  
blood loss

Side port and  
3-way stopcock  
for convenient flushing

Braided sheath  
for optimal flexibility  
and kink resistance




Radiopaque  
marker band  
for easy imaging



**Available in  
30 cm, 65 cm,  
and 85 cm  
lengths**

## FOR YOUR MOST COMPLEX CASES

The D'VILL™ is made of an innovative braided material for enhanced flexibility and kink resistance. Available in 12F and 14F sizes, the D'VILL™ is compatible with CP stents® and BIB® catheters.



### Indications for use:

Recommended for introduction of balloons, catheters, and other diagnostic and interventional devices. Refer to the D'VILL™ Instructions for Use for relevant warnings, precautions, complications, and contraindications.

# D'VILL™ Introducer

Long-length, kink-resistant  
hemostasis introducer



## For Introduction of Large Diameter Interventional Devices

Compatible with CP Stent® and BIB® catheters

D'VILL™ Specifications			
Catalog No.	Introducer Size (FR)	Usable Length (CM)	Guidewire (inches)
DV1230	12	30	0.035
DV1430	14	30	0.035
DV1265	12	65	0.035
DV1465	14	65	0.035
DV1285	12	85	0.035
DV1485	14	85	0.035

- Each D'VILL™ set includes (1) one sheath and (1) one dilator

The D'VILL™ is available worldwide. Please contact us or your local distributor to place an order.

*This device is subject to individual country regulations regarding the importation and sale of this product*



**WORLD LEADER IN  
PEDIATRIC CARDIOLOGY**

NuMED | 2880 Main Street  
Hopkinton, NY 12965 USA  
Tel: 315.328.4491 | Fax: 315.328.4941  
[www.numedforchildren.com](http://www.numedforchildren.com)

Rev\_2021-06-01\_BVM

**BVM Medical Ltd.**  
UK Distributor  
[info@bvmmedical.com](mailto:info@bvmmedical.com)  
BVM House, Trinity Lane  
Hinckley, Leicestershire  
LE10 0BL United Kingdom  
Tel: **01455-614555**  
Fax: **01455-614546**

**BVM**  
— 24 hours a day —

**G. van Wageningen B.V.**  
European Representative  
[info@heartmedical.nl](mailto:info@heartmedical.nl)  
Hallenweg 40  
5683 CT Best  
The Netherlands  
Tel: **+31-499-377388**  
Fax: **+31-499-377456**

**CE 1639**