


01



EGIS Flower Stent

- Effective drainage of bile juice by side grooves
- Prevent cholecystitis and pancreatitis
- e-PTFE fully covered and retrieval stent

Upgraded Special Features



- Innovative design with 5-6 grooves, facilitating bile drainage



- Fully covered stent with e-PTFE, designed for use in malignant & benign strictures



- Sufficient radial force to push out the stricture



- Low axial force, ensuring excellent flexibility & conformability



- For easy removal, a triangle-shaped lasso is attached



- To obtain anti-migration function, flared ends and grooves are made in the entire body

Delivery System 8Fr, 50/180cm

Ordering NO. First two** : Diameter(mm) / Last two** : Length(cm)

Ordering No.			Size	
Type	OTW	TTS	Diameter	Length
Bare	BFBB****0805	BFBB****0818	8, 10	4, 5, 6, 7, 8, 9, 10, 12
Partial	BFPB****0805	BFPB****0818		
Fully Covered	BFCB****0805	BFCB****0818		

Stent Type	OTW(Over The Wire)					TTS(Through The Scope)				
	Ordering NO.	Stent		Delivery System		Ordering No.	Stent		Delivery System	
		Diameter(mm)	Length(cm)	Profile(fr.)	Length(mm)		Diameter(mm)	Length(cm)	Profile(fr.)	Length(mm)
Flower Fully Covered	BFCB08040805	8	4	8	500	BFCB08040818	8	4	8	1800
	BFCB08050805		5			BFCB08050818		5		
	BFCB08060805		6			BFCB08060818		6		
	BFCB08070805		7			BFCB08070818		7		
	BFCB08080805		8			BFCB08080818		8		
	BFCB08090805		9			BFCB08090818		9		
	BFCB08100805		10			BFCB08100818		10		
	BFCB08120805		12			BFCB08120818		12		
	BFCB10040805	10	4	8	500	BFCB10040818	10	4	8	1800
	BFCB10050805		5			BFCB10050818		5		
	BFCB10060805		6			BFCB10060818		6		
	BFCB10070805		7			BFCB10070818		7		
	BFCB10080805		8			BFCB10080818		8		
	BFCB10090805		9			BFCB10090818		9		
	BFCB10100805		10			BFCB10100818		10		
	BFCB10120805		12			BFCB10120818		12		
Flower Partial Covered	BFPB08040805	8	4	8	500	BFPB08040818	8	4	8	1800
	BFPB08050805		5			BFPB08050818		5		
	BFPB08060805		6			BFPB08060818		6		
	BFPB08070805		7			BFPB08070818		7		
	BFPB08080805		8			BFPB08080818		8		
	BFPB08090805		9			BFPB08090818		9		
	BFPB08100805		10			BFPB08100818		10		
	BFPB08120805		12			BFPB08120818		12		
	BFPB10040805	10	4	8	500	BFPB10040818	10	4	8	1800
	BFPB10050805		5			BFPB10050818		5		
	BFPB10060805		6			BFPB10060818		6		
	BFPB10070805		7			BFPB10070818		7		
	BFPB10080805		8			BFPB10080818		8		
	BFPB10090805		9			BFPB10090818		9		
	BFPB10100805		10			BFPB10100818		10		
	BFPB10120805		12			BFPB10120818		12		
Flower Bare	BFBB08040805	8	4	8	500	BFBB08040818	8	4	8	1800
	BFBB08050805		5			BFBB08050818		5		
	BFBB08060805		6			BFBB08060818		6		
	BFBB08070805		7			BFBB08070818		7		
	BFBB08080805		8			BFBB08080818		8		
	BFBB08090805		9			BFBB08090818		9		
	BFBB08100805		10			BFBB08100818		10		
	BFBB08120805		12			BFBB08120818		12		
	BFBB10040805	10	4	8	500	BFBB10040818	10	4	8	1800
	BFBB10050805		5			BFBB10050818		5		
	BFBB10060805		6			BFBB10060818		6		
	BFBB10070805		7			BFBB10070818		7		
	BFBB10080805		8			BFBB10080818		8		
	BFBB10090805		9			BFBB10090818		9		
	BFBB10100805		10			BFBB10100818		10		
	BFBB10120805		12			BFBB10120818		12		