

TECHNICAL DETAILS


INTRODUCER SHEATH		DILATOR	
Material – outer	Stainless steel reinforced polymer	Tip	Smooth taper, radiopaque
Material – inner	PTFE liner	Shape	Straight
Sheath coating	Hydrophobic, distal 15 cm		
Sheath tip	Smooth taper with radiopaque platinum iridium marker band		
French size	4, 5, 6, 7, 8 F		
Length	15, 45, 60, 90 cm		
Sidearm/hemostatic valve	Removable, with color coded 3-way-stopcock		

ORDERING INFORMATION

CATALOGUE NUMBER	DIAMETER (FRENCH)	LENGTH (CM)	INNER DIAMETER MINIMUM (MM/INCH)	DILATOR GUIDEWIRE COMPATIBILITY MAX.	UNITS/ CASE
FG-04015	4 F	15 cm	1.60 MM / 0.063"	0.035" & 0.018"	5
FG-04045	4 F	45 cm	1.60 MM / 0.063"	0.035"	5
FG-04060	4 F	60 cm	1.60 MM / 0.063"	0.035"	5
FG-04090	4 F	90 cm	1.60 MM / 0.063"	0.035"	5
FG-05015	5 F	15 cm	1.92 MM / 0.075"	0.035" & 0.018"	5
FG-05045	5 F	45 cm	1.92 MM / 0.075"	0.035"	5
FG-05060	5 F	60 cm	1.92 MM / 0.075"	0.035"	5
FG-05090	5 F	90 cm	1.92 MM / 0.075"	0.035"	5
FG-06015	6 F	15 cm	2.25 MM / 0.088"	0.035" & 0.018"	5
FG-06045	6 F	45 cm	2.25 MM / 0.088"	0.035"	5
FG-06060	6 F	60 cm	2.25 MM / 0.088"	0.035"	5
FG-06090	6 F	90 cm	2.25 MM / 0.088"	0.035"	5
FG-07045	7 F	45 cm	2.55 MM / 0.100"	0.038"	5
FG-07060	7 F	60 cm	2.55 MM / 0.100"	0.038"	5
FG-07090	7 F	90 cm	2.55 MM / 0.100"	0.038"	5
FG-08045	8 F	45 cm	2.88 MM / 0.113"	0.038"	5
FG-08060	8 F	60 cm	2.88 MM / 0.113"	0.038"	5
FG-08090	8 F	90 cm	2.88 MM / 0.113"	0.038"	5

15 cm versions contain two dilators.
Additional sizes to become available.

Manufactured by:
Contract Medical International GmbH
Lauensteiner Str. 37
01277 Dresden
Germany



KINK-RESISTANT
INTRODUCER SHEATH
Supporting complex access

SUPPORTING COMPLEX ACCESS

From easy insertion to smooth navigation to superb visibility:
DuraSheath is designed to combine optimal performance with excellent pushability and trackability. While maintaining the inner diameter, DuraSheath allows for smooth device passage and has an exceptional kink-radius and non-ovaling effect. The detachable hemostatic valve provides flexibility and versatility, giving you the ability to focus on the intervention.

