



Product Catalog

Mongoose[®] Pediatric Angiography Catheters
Super Sheath[®] Pediatric Introducers
Spring Guide Wires

CE MARKED FOR PEDIATRIC USE

BVM
24 hours a day

01455 614555
www.bvmmedical.com
info@bvmmedical.com

Mongoose®

Pediatric Angiography Catheters

CE MARKED For Pediatric Use

Product Code	Size	Tip Shape	Length	Maximum Wire	Maximum Flow Rate*	
A-3PIG-40/1001	3.3F	Pigtail	40 cm	.030"	13 mL/sec	
A-3PIG-60/1002	3.3F	Pigtail	60 cm	.030"	11 mL/sec	
A-3PIG-80/1003	3.3F	Pigtail	80 cm	.030"	9 mL/sec	
A-4PIG-60/1004	4F	Pigtail	60 cm	.035"	22 mL/sec	
A-4PIG-80/1005	4F	Pigtail	80 cm	.035"	19 mL/sec	
A-3JB1-60/1008	3.3F	JB1	60 cm	.030"	11 mL/sec	
A-3JB1-80/1009	3.3F	JB1	80 cm	.030"	9 mL/sec	
A-3JB1-60/1044 ¹	3.3F	JB1	60 cm	.030"	11 mL/sec	
A-3JB1-80/1045 ¹	3.3F	JB1	80 cm	.030"	9 mL/sec	
A-4JB1-60/1010	4F	JB1	60 cm	.038"	22 mL/sec	
A-4JB1-80/1011	4F	JB1	80 cm	.038"	19 mL/sec	
A-3JR1-60/1014	3.3F	JR1	60 cm	.030"	11 mL/sec	
A-3JR1-80/1015	3.3F	JR1	80 cm	.030"	9 mL/sec	
A-4JR1-60/1016	4F	JR1	60 cm	.038"	22 mL/sec	
A-4JR1-80/1017	4F	JR1	80 cm	.038"	19 mL/sec	
A-3JR2-60/1046	3.3F	JR2	60 cm	.030"	11 mL/sec	
A-3JR2-80/1047	3.3F	JR2	80 cm	.030"	9 mL/sec	
A-4JR2-60/1048	4F	JR2	60 cm	.038"	22 mL/sec	
A-4JR2-80/1049	4F	JR2	80 cm	.038"	19 mL/sec	

¹ No side holes, designed for placing coils

*@ 1200 PSI

Mongoose®

Pediatric Angiography Catheters

CE MARKED For Pediatric Use


Product Code	Size	Tip Shape	Length	Maximum Wire	Maximum Flow Rate*	
A-3JL1-60/1020	3.3F	JL1	60 cm	.030"	11 mL/sec	
A-3JL1-80/1021	3.3F	JL1	80 cm	.030"	9 mL/sec	
A-4JL1-60/1022	4F	JL1	60 cm	.038"	22 mL/sec	
A-4JL1-80/1023	4F	JL1	80 cm	.038"	19 mL/sec	
A-3JLA-60/1026	3.3F	JL1.5	60 cm	.030"	11 mL/sec	
A-3JLA-80/1027	3.3F	JL1.5	80 cm	.030"	9 mL/sec	
A-4JLA-60/1028	4F	JL1.5	60 cm	.038"	22 mL/sec	
A-4JLA-80/1029	4F	JL1.5	80 cm	.038"	19 mL/sec	
A-3MP-60/1032	3.3F	MP	60 cm	.030"	11 mL/sec	
A-3MP-80/1033	3.3F	MP	80 cm	.030"	9 mL/sec	
A-4MP-60/1034	4F	MP	60 cm	.038"	22 mL/sec	
A-4MP-80/1035	4F	MP	80 cm	.038"	19 mL/sec	
A-3COB-60/1038	3.3F	Cobra	60 cm	.030"	11 mL/sec	
A-3COB-80/1039	3.3F	Cobra	80 cm	.030"	9 mL/sec	
A-4COB-60/1040	4F	Cobra	60 cm	.038"	22 mL/sec	
A-4COB-80/1041	4F	Cobra	80 cm	.038"	19 mL/sec	

*@ 1200 PSI

Super Sheath[®]

Pediatric Introducers

FDA Cleared For Pediatric Use



Product Code	Size	Dilator Taper	Length	Wire (kit)
S-3N5/1005	3.3F	.018"	5 cm	n/a
S-3N7/1006	3.3F	.018"	7 cm	n/a
S-3W5/1007	3.3F	.018"	5 cm	.018"
S-3W7/1008	3.3F	.018"	7 cm	.018"

Spring Wires

Product Code	Size	Length	Floppy Tip
SW30-140	.030"	140 cm	2 cm
SW30-140J	.030"	140 cm	5 mm J

BVM

24 hours a day

01455 614555
www.bvmmedical.com
info@bvmmedical.com